

Your teacher may watch to see if you can...

- safely and correctly use apparatus.

Aim

To prepare a sample of pure, dry, hydrated copper sulfate crystals starting from copper oxide.

Apparatus

- | | |
|--|-----------------------------|
| • eye protection | • spatula |
| • 100 cm ³ conical flask | • stirring rod |
| • 100 cm ³ beaker | • filter funnel |
| • Bunsen burner | • filter paper |
| • gauze and tripod | • tongs |
| • heat mat | • water bath (set at 50 °C) |
| • Petri dish or watch glass | • dilute sulfuric acid |
| • 100 cm ³ measuring cylinder | • copper(II) oxide |
| • evaporating basin | |

Safety

Wear eye protection at all times.

Method

- A** Pour about 20 cm³ of dilute sulfuric acid into a conical flask.
- B** Place the conical flask into a water bath at 50 °C and heat for 3–4 minutes to allow the acid to heat up.
- C** Use the spatula to add a little copper oxide to the acid and stir or swirl the contents of the flask.
- D** Keep repeating step **C** until the black powder does not disappear after stirring. (This makes sure the copper oxide is in excess.)
- E** Return the mixture to the water bath for a few minutes (to make sure there is no more acid left).
- F** Filter the mixture into a beaker and pour into an evaporating basin.
- G** Place the evaporating basin on top of a beaker half full of water. Heat the beaker, evaporating basin and contents using a Bunsen burner on a blue flame.
- H** Heat until about half of the water has evaporated. Then allow the evaporating basin to cool.
- I** When cool, transfer the solution to a Petri dish or watch glass and leave for a few days to allow the water to evaporate.
- J** Observe the shape and colour of the copper sulfate crystals formed.



Step F



Step H